



Artificial Intelligence Algorithm Detecting Status Epilepticus and Measuring Seizure Burden

Baharan Kamousi PhD, Archit Gupta PhD, Suganya Karunakaran PhD, Ali Marjaninejad PhD, Ray Woo PhD, Josef Parvizi MD, PhD

INTRODUCTION

- Artificial intelligence (AI) algorithms are increasingly integrated into medical practice – improving efficiency, quality of life, and diagnostic decision making of clinical decision makers.
- The FDA-cleared **Ceribell Clarity Algorithm** automatically and continuously monitors the EEG, measuring the *seizure burden* every 10 seconds (i.e., burden of seizure activity in the last 5 minutes)
- When suspected status epilepticus is detected, it provides visual and auditory alerts at the bedside and remotely.
- With added data, AI algorithms have the potential to continuously improve.
- The current study was a comparative study to measure the accuracy of the latest version (v6.0) of the algorithm compared to the earlier version of the same algorithm (v2.0)

METHODS

- We retrospectively selected 666 Ceribell EEGs from 11 centers.
- Each EEG was categorized using the majority consensus of at least two expert epileptologists.

RESULTS


- The majority Epileptologist consensus labeled the EEGs as following:

Slow/Non-epileptiform	Highly Epileptiform	Seizure	Status Epilepticus
501	128	17	20


- The Clarity AI output was as follows:

Clarity v6.0 Seizure Burden	Human Epileptologist Rating			
	Slow/Non-epileptiform	Highly Epileptiform	Seizure	Status Epilepticus
0%	392	54	4	0
>0%, <90%	108	59	12	1
≥90%	1	15	1	19
Total	501	128	17	20

Clarity analyzes the EEG every 10 seconds to output seizure burden



Seizure burden display and alert on the Ceribell device.



Clarity v6.0 Algorithm Accuracy

95%	Sensitivity Accurately detecting status epilepticus
97%	Specificity Accurately identifying non-status cases as not status epilepticus
99.8%	Negative Predictive Value Accurately ruling out status epilepticus

- In the single missed case of status epilepticus, Clarity still indicated a high degree of seizure activity.
- Clarity had positive predictive value (PPV) of **53%** for status, due to overcalling 17 of 666 recordings. All but one of the 17 overcalled recordings had been labeled by majority consensus of epileptologists as highly epileptiform or seizure.
- There was a large improvement in the positive predictive value of the ≥90% seizure burden to alert to status epilepticus (32% to 53%) as the algorithm alerted to half as many non-SE cases compared to the earlier version.
- The negative predictive value (NPV) for 90% seizure burden ruling out only the presence of status epilepticus remained unchanged at 99.8% while the NPV for 0% seizure burden ruling out the presence of status epilepticus or any seizure improved from 98.4% to 99.1%

Comparative data for Seizure Burden ≥90%:

Clarity Algorithm	Sensitivity	Specificity	PPV	NPV
v2.0	95%	94%	32%	99.8%
v6.0	95%	97%*	53%	99.8%

*p<0.001

Conclusions

- This study demonstrates the ability of machine learning algorithms to improve with time and additional cases for training.
- The latest algorithm led to an improvement in specificity without sacrificing high sensitivity to SE.
- The high negative predictive value of the algorithm at 0% threshold suggests that cases of status epilepticus can be ruled out relatively accurately in a large proportion of cases within minutes of EEG recordings and thus can help prevent unnecessary or aggressive over-treatment in critical care settings – as shown in recent clinical studies.